

## IMPORTANT: Please read carefully and keep this information for future use.

### Boston SIMPLUS™ Multi-Action Solution

#### For Cleaning, Removing Protein, Rinsing, Disinfecting and Conditioning Fluoro Silicone Acrylate and Silicone Acrylate Rigid Gas Permeable Contact Lenses

**Description:** A sterile, aqueous, buffered solution that contains poloxamine, hydroxyalkylphosphonate, boric acid, sodium borate, sodium chloride, hydroxypropylmethyl cellulose, Glucam and preserved with chlorhexidine gluconate (0.003%), polyaminopropyl biguanide (0.0005%).

**Actions:** This solution removes protein and lipids from the lens surface and also helps keep stubborn deposits away from lenses. This solution is also designed for rinsing lenses, which eliminates loosened surface deposits from the lens. A multi-polymer cushioning system conditions the lens to provide cushioning upon insertion. A patented humectant maintains lens wettability. This formula is designed to provide a soothing effect to your eyes. It also contains a preservative system proven effective in destroying harmful organisms. This solution is effective both in cleaning and disinfecting, yet is compatible with the eye upon lens insertion.

**Indications:** For cleaning, removing protein, rinsing, disinfecting, conditioning, storing and cushioning fluoro silicone acrylate and silicone acrylate rigid gas permeable contact lenses.

**Contraindications:** Do not use this solution if you are allergic to any ingredient in this product. Not For Use With Soft (Hydrophilic) Contact Lenses.

**Warnings:** To avoid contamination, never touch dropper tip of the container with your hands or to any surface. Keep bottle tightly closed when not in use. Contact lenses and solutions, if not handled properly, can become contaminated by common bacteria known to cause ocular infections. Serious damage to the eye and loss of vision may result from using contaminated solutions or from wearing contaminated or damaged contact lenses.

**PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO YOUR EYES. It is essential that you follow your eye care professional's directions and all labeling instructions for proper use of your lenses and lens care products.**

**EYE PROBLEMS, INCLUDING CORNEAL ULCERS CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISION. THEREFORE, IF YOU EXPERIENCE EYE DISCOMFORT, EXCESSIVE TEARING, VISION CHANGES OR REDNESS OF THE EYES, IMMEDIATELY REMOVE YOUR LENSES AND PROMPTLY CONTACT YOUR EYE CARE PRACTITIONER.**

**Precautions:** To ensure proper hygiene and handling of contact lenses and contact lens solutions, always follow label directions carefully.

- Never reuse this solution.
- Always use fresh solution for soaking and storing lenses.
- Always wash your hands before handling your lenses.
- Store solutions at room temperature. Avoid freezing.
- After inserting your lenses, always thoroughly clean the interior of the lens case with hot water and air dry after each use.
- Use before expiration date on the container.
- ALWAYS RINSE LENSES WITH BOSTON SIMPLUS MULTI-ACTION SOLUTION AFTER CLEANING TO ELIMINATE LOOSENED SURFACE DEPOSITS.
- KEEP OUT OF THE REACH OF CHILDREN.





**Adverse Effects: The following problems may occur:**

- Eyes stinging, burning, itching (irritation)
- Excessive watering (tearing) of the eyes
- Unusual eye secretions
- Redness of the eyes
- Reduced sharpness of vision (visual quality)
- Blurred vision
- Sensitivity to light (photophobia)
- Dry eyes

If you notice any of the above conditions, immediately remove and examine the lenses. If the problem stops and the lenses appear to be undamaged, thoroughly clean, rinse and disinfect the lenses and reinsert. If the problem continues or if a lens appears to be damaged, DO NOT REINSERT. Contact your eye care professional immediately. If any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and treatment to avoid serious eye damage.

### DIRECTIONS FOR USE

To ensure proper disinfecting, all steps listed below must be followed:

-  Wash your hands with mild soap. (Caution - pump soaps may contain oil-based suspension agents.)
-  Place lenses in empty lens case and fill with fresh Boston SIMPLUS Multi-Action Solution. **Soak lenses for at least four (4) hours (or overnight) before wearing.**
-  Wash hands with mild soap before cleaning (rubbing) lenses. After soaking, remove lenses from lens case and rub both sides of the lenses carefully with four (4) drops of Boston SIMPLUS Multi-Action Solution in the palm of your hand for twenty (20) seconds. **NO SEPARATE DAILY CLEANER IS REQUIRED.**
-  Rinse for approximately five (5) seconds with a steady stream of Boston SIMPLUS Multi-Action Solution to eliminate loosened surface deposits and insert lenses.
5. Thoroughly clean interior of lens case with hot water and air dry after each use.
6. Discard solution ninety (90) days after opening. Record date opened in space provided on bottle label. Note: When used daily, this solution will be depleted before ninety (90) days.

- For lens rewetting and lubrication, use Boston® Rewetting Drops
- Lens cases should be replaced frequently to avoid potential microbial contamination and ocular irritation.

**How Supplied:** Boston SIMPLUS Multi-Action Solution is supplied sterile in 3 fluid ounce (90 mL) and 3.5 fluid ounce (105 mL) bottles. All retail Boston® Solutions are packed in unit cartons. Bottles are sealed with an imprinted neckband and marked with a lot number and expiration date. **Do not use if imprinted neckband on bottle is missing or broken.**

© Polymer Technology, a Bausch & Lomb company  
Made in USA 4/03

Manufactured for:  
Polymer Technology,  
a Bausch & Lomb company,  
Rochester, NY 14609

Questions or Comments:  
**1-800-333-4730**  
[www.bausch.com](http://www.bausch.com)